



Huvitz
Pacing Progress toward People



HCP-7000 CHART PROJECTOR

- Semi-permanent lifetime LED light source
- Brighter and clearer than conventional bulb light
- Attractive Streamlined Shape and Blue Indicator
- RC-R Type (Russian Chart) added

SPECIFICATIONS

Chart	41 charts 34 masks, Red/Green & Polarization Filters
Projection Distance	2.5 ~ 8m
Chart Rotation Speed	Average 0.15 sec
Projection Magnification	30x at 5m
Power Saving	Automatic Switch off (10 min)
Program	2 programs with a maximum of 30 charts each
Tilt Angle	15°
Power Supply	100-120V 50Hz : 0.6A, 200-240V 60Hz : 0.3A
Lamp	LED 4W
Dimension	270(W) x 182(D) x 230(H)mm
Weight	3.44kg

Designs and details can be changed without prior notice for the purposes of improvement.



9990 NW 14th Street
Unit # 105
Doral, Florida 33172
Tel: 305-722-7321
Fax: 786-664-3347
info@usophthalmic.com
www.usophthalmic.com

Distributed by